

Folic acid Monoclonal Antibody



MACO0590

Product Information

Size:

50ug

Reactivity:

Source:

Mouse

Isotype:

IgG2b

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Folic acid

Uniprot

Synonyms:

Acidum folicum, aminomethopterin, folacin, folic acid, pteroylglutamic acid, pteroylmonoglutamic acid, vitamin B9, vitamin L, vitamin M, PteGlu, PTGA, FA, VM

Immunogen:

Folic acid-BSA conjugate

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A