

Sulfametoxydiazine Monoclonal Antibody



MACO0591

Product Information

Size:

50ug

Reactivity:

Source:

Mouse

Isotype:

IgG1

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Sulfametoxydiazine

Uniprot

Synonyms:

SMD, aristebon, sulfametoxydiazine, sulfameter, bayrena, durenat, 5-Methoxysulfadiazine

Immunogen:

Sulfametoxydiazine-BSA conjugate

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A