

Sulfadimidine Monoclonal Antibody



MACO0592

Product Information

Size:

50ug

Reactivity:

Source:

Mouse

Isotype:

IgG2a

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Sulfadimidine

Uniprot

Synonyms:

SM2, S-Mez, Sulfamethazine, 4,6-Dimethylsulfadiazine, 4-Amino-N-(4,6-dimethyl-2-pyrimidinyl)benzenesulf-onamide

Immunogen:

Sulfadimidine-BSA conjugate

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A