

estradiol Monoclonal Antibody



MACO0595

Product Information

Size:

50ug

Protein Background:**Reactivity:****Gene ID:**

estradiol

Source:**Uniprot**

Mouse

Isotype:**Synonyms:**

IgG2b

Beta-estradiol, 17beta-estradiol, estrace, oestradiol, dihydrofolliculin, progynon, dihydroxyestrin, dihydrotheelin, diogynets

Applications:**Immunogen:**

ELISA

estradiol acid-BSA conjugate

Recommended dilutions:**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A