

# Cortisol Monoclonal Antibody



MACO0596

---

## Product Information

**Size:**

50ug

**Reactivity:**

**Source:**

Mouse

**Isotype:**

IgG1

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

Cortisol

**Uniprot**

**Synonyms:**

Hydrocortisone

**Immunogen:**

Cortisol-BSA conjugate

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A