

abscisic acid Monoclonal Antibody



MACO0601

Product Information

Size:

50ug

Protein Background:**Reactivity:****Gene ID:**

abscisic acid

Source:**Uniprot**

Mouse

Isotype:**Synonyms:**

IgG1

ABA, abscissic acid, dormin

Applications:**Immunogen:**

ELISA

abscisic acid-BSA conjugate

Recommended dilutions:**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A