

Testosterone Monoclonal Antibody



MACO0602

Product Information

Size:

50ug

Reactivity:

Source:

Mouse

Isotype:

IgG2b

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Testosterone

Uniprot

Synonyms:

Ts, andrusol, sustanon, teslen, testicular hormone, virosterone

Immunogen:

Testosterone-BSA conjugate

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A