

Cleaved-C1S (R437) Antibody



PACO00008

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

C1S

Uniprot

P09871

Synonyms:

C1S; Complement C1s subcomponent; C1 esterase; Complement component 1 subcomponent s

Immunogen:

Synthesized peptide derived from the Internal region of human C1s HC.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using Cleaved-C1s HC (R437) Polyclonal Antibody.