

Cleaved-CTSZ (L62) Antibody



PACO00025

Product Information

Size:

50ug

Reactivity:

Human, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

CTSZ

Uniprot

Q9UBR2

Synonyms:

CTSZ; Cathepsin Z; Cathepsin P; Cathepsin X


Immunogen:

Synthesized peptide derived from the Internal region of human Cathepsin Z.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COS7 cells using Cleaved-Cathepsin Z (L62) Polyclonal Antibody.