

# Cleaved-COL3A1 (G1221) Antibody



PACO00026


---

## Product Information

<b>Size:</b>	<b>Protein Background:</b>
50ug	
<b>Reactivity:</b>	<b>Gene ID:</b>
Human	COL3A1
<b>Source:</b>	<b>Uniprot</b>
Rabbit	P02461
<b>Isotype:</b>	<b>Synonyms:</b>
IgG	COL3A1; Collagen alpha-1(III) chain
<b>Applications:</b>	<b>Immunogen:</b>
ELISA, WB	Synthesized peptide derived from the C-terminal region of human COL3A1.
<b>Recommended dilutions:</b>	<b>Storage:</b>
WB:1:500-1:2000	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of A549 cells using Cleaved-COL3A1 (G1221) Polyclonal Antibody.