

Cleaved-LMNA (D230) Antibody



PACO00032

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

LMNA

Uniprot

P02545

Synonyms:

LMNA; LMN1; Prelamin-A/C


Immunogen:

Synthesized peptide derived from the Internal region of human Lamin A.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using Cleaved-Lamin A (D230) Polyclonal Antibody.