

Cleaved-MMP23A (Y79) Antibody




PACO00036

Product Information

Size: 50ug	Protein Background:
Reactivity: Human, Mouse, Rat	Gene ID: MMP23A
Source: Rabbit	Uniprot O75900
Isotype: IgG	Synonyms: MMP23A; MMP21; MMP23B; MMP21; MMP22; Matrix metalloproteinase-23; MMP-23; Femalysin; MIFR-1; Matrix metalloproteinase-21; MMP-21; Matrix metalloproteinase-22; MMP-22
Applications: ELISA, WB	Immunogen: Synthesized peptide derived from the Internal region of human MMP-23.
Recommended dilutions: WB:1:500-1:2000	Storage: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293 cells using Cleaved-MMP-23 (Y79) Polyclonal Antibody.