

Cleaved-MMP27 (Y99) Antibody




PACO00037

Product Information

Size: 50ug	Protein Background:
Reactivity: Human, Monkey	Gene ID: MMP27
Source: Rabbit	Uniprot Q9H306
Isotype: IgG	Synonyms: MMP27; Matrix metalloproteinase-27; MMP-27
Applications: ELISA, WB	Immunogen: Synthesized peptide derived from the Internal region of human MMP-27.
Recommended dilutions: WB:1:500-1:2000	Storage: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COS7 cells using Cleaved-MMP-27 (Y99) Polyclonal Antibody.