

Cleaved-F10 (A41) Antibody



PACO00044

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

F10

Uniprot

P00742

Synonyms:

F10; Coagulation factor X; Stuart factor; Stuart-Prower factor

Immunogen:

Synthesized peptide derived from the N-terminal region of human Factor X/ Factor X LC.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of A549 cells using Cleaved-Factor X/ Factor X LC (A41) Polyclonal Antibody.