

Acetyl-ARNTL (K538) Antibody



PACO00069

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

ARNTL

Uniprot

O00327

Synonyms:

ARNTL; BHLHE5; BMAL1; MOP3; PASD3; Aryl hydrocarbon receptor nuclear translocator-like protein 1; Basic-helix-loop-helix-PAS protein MOP3


Immunogen:

Synthesized peptide derived from human BMAL1 around the acetylation site of K538.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293 cells using Acetyl-BMAL1 (K538) Polyclonal Antibody.