

Acetyl-TUBA1A/TUBA1B/TUBA1C (K112)
Antibody
PACO00072



Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

TUBA1A/TUBA1B/TUBA1C

Uniprot

Q71U36/P68363/Q9BQE3

Synonyms:

TUBA1A; TUBA3; Tubulin alpha-1A chain; Alpha-tubulin 3; Tubulin B-alpha-1; Tubulin alpha-3 chain; TUBA1B; Tubulin alpha-1B chain; Alpha-tubulin ubiquitous


Immunogen:

Synthesized peptide derived from human Tubulin alpha around the acetylation site of K112.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HepG2 cells using Acetyl-Tubulin alpha (K112) Polyclonal Antibody.