

Myc-Tag Antibody



PACO00078

Product Information

Size:

50ug

Protein Background:**Reactivity:****Gene ID:****Source:**

Rabbit

Uniprot**Isotype:**

IgG

Synonyms:**Applications:**

WB


Immunogen:

Synthetic peptide.

Recommended dilutions:**Storage:**

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



1ug C-Myc fusion protein+ Primary antibody dilution at 1) 1:5000 2) 1:10000.