

TUBA1A/TUBA1B Antibody



PACO00086

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

WB

Recommended dilutions:

Protein Background:

Gene ID:

TUBA1A/TUBA1B

Uniprot

Q71U36/P68363

Synonyms:

Tubulin alpha-1A chain


Immunogen:

Recombinant Protein.

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of 1) Mouse skeletal muscle tissue, 2) Mouse kidney tissue, 3) Mouse Heart Tissue, 4) Mouse brain tissue, 5) Mouse Liver tissue, diluted at 1:2000.