

MECP2 Antibody



PACO00088

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

WB

Recommended dilutions:

Protein Background:

Gene ID:

MECP2

Uniprot

P51608

Synonyms:

Methyl-CpG-binding protein 2


Immunogen:

Synthetic peptide.

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of HeLa, diluted at 1) 1:2000 2) 1:5000.