

CREB1 Antibody



PACO00099

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

WB

Recommended dilutions:**Protein Background:****Gene ID:**

CREB1

Uniprot

P16220

Synonyms:

CREB1; Cyclic AMP-responsive element-binding protein 1; CREB-1; cAMP-responsive element-binding protein 1


Immunogen:

Recombinant Protein.

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of HeLa, Rabbit pAb diluted at 1:2000.