

BSA Antibody



PACO00112

Product Information

Size:

50ug

Reactivity:

Cow

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:2000-1:5000

Protein Background:

Gene ID:

Uniprot

Synonyms:


Immunogen:

Recombinant Protein.

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of Bovine Serum Albumin, Rabbit pAb diluted at
1) 1:2000 2) 1:5000.