

# Phospho-Histone H2A (Ser129) Antibody



PACO00136

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

WB

**Recommended dilutions:****Protein Background:****Gene ID:**

HIST1H2AG/HIST2H2AA3/HIST3H2A

**Uniprot**

P0C0S8/Q6FI13/Q7L7L0

**Synonyms:**

HIST1H2AG; H2AFP; HIST1H2AI; H2AFC; HIST1H2AK; H2AFD; HIST1H2AL; H2AFI; HIST1H2AM; H2AFN; Histone H2A type 1; H2A.1; Histone H2A/p; HIST2H2AA3; H2AFO; HIST2H2AA; HIST2H2AA4; Histone H2A type 2-A; Histone H2A.2; Histone H2A/o; HIST3H2A; Histone H2A type 3

**Immunogen:**


Synthetic peptide.

**Storage:**

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Product Images

---



Western blot analysis of extracts from HeLa cells, untreated (-) or treated, 1:5000.