

# Acetyl-Histone H2A. Z (Lys7) Antibody



PACO00158

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:1000

**Protein Background:****Gene ID:**

H2AFZ

**Uniprot**

P0C0S5

**Synonyms:**

H2AFZ; H2AZ; Histone H2A. Z; H2A/z

**Immunogen:**


Synthetic peptide.

**Storage:**

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Product Images

---



Western blot analysis of extracts from HeLa cells, untreated (-) or treated, 1:5000.