

# Phospho-Histone H2A. X (Thr120) Antibody



PACO00169

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

H2AFX

**Uniprot**

P16104

**Synonyms:**

H2AFX; H2AX; Histone H2A. x; H2a/x

**Immunogen:**


Synthetic peptide.

**Storage:**

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Product Images

---



Western blot analysis of extracts from Hela cells, 1:2000.