

Acetyl-Histone H2A. X (Lys5) Antibody



PACO00173

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

WB

Recommended dilutions:

Protein Background:

Gene ID:

H2AFX

Uniprot

P16104

Synonyms:

H2AFX; H2AX; Histone H2A. x; H2a/x


Immunogen:

Synthetic peptide.

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of extracts from Hela cells, untreated (-) or treated, 1:5000.