

Phospho-PPP1R1B (T75) Antibody



PACO00199

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

PPP1R1B

Uniprot

Q9UD71

Synonyms:

PPP1R1B; DARPP32; Protein phosphatase 1 regulatory subunit 1B; DARPP-32; Dopamine- and cAMP-regulated neuronal phosphoprotein


Immunogen:

Synthesized peptide derived from human DARPP-32 around the phosphorylation site of T75.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using Phospho-DARPP-32 (T75) Polyclonal Antibody.