

# Phospho-EPHA2 (Y588) Antibody



PACO00206

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:****Gene ID:**

EPHA2

**Uniprot**

P29317

**Synonyms:**

EPHA2; ECK; Ephrin type-A receptor 2; Epithelial cell kinase; Tyrosine-protein kinase receptor ECK

**Immunogen:**

Synthesized peptide derived from human EphA2 around the phosphorylation site of Y588.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of HepG2 cells using Phospho-EphA2 (Y588) Polyclonal Antibody.

