

Phospho-ESR1 (S104) Antibody



PACO00208

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

ESR1

Uniprot

P03372

Synonyms:

ESR1; ESR; NR3A1; Estrogen receptor; ER; ER-alpha; Estradiol receptor; Nuclear receptor subfamily 3 group A member 1

Immunogen:

Synthesized peptide derived from human ERalpha around the phosphorylation site of S104.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of 293 cells using Phospho-ERalpha (S104) Polyclonal Antibody.

