

Phospho-FOXO4 (S197) Antibody



PACO00211

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

FOXO4

Uniprot

P98177

Synonyms:

FOXO4; AFX; AFX1; MLLT7; Forkhead box protein O4; Fork head domain transcription factor AFX1


Immunogen:

Synthesized peptide derived from human FoxO4 around the phosphorylation site of S197.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using Phospho-FoxO4 (S197) Polyclonal Antibody.