

Phospho-CHUK/IKBKB (S180/181) Antibody



PACO00221

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

CHUK/IKBKB

Uniprot

O15111/O14920

Synonyms:

CHUK; IKKA; TCF16; Inhibitor of nuclear factor kappa-B kinase subunit alpha; I-kappa-B kinase alpha; IKK-A; IKK-alpha; Ikbka; Ikbkba; Ikbkba kinase; Conserved helix-loop-helix ubiquitous kinase; I-kappa-B kinase 1; IKK1; Nuclear factor NF-kappa-B


Immunogen:

Synthesized peptide derived from human IKKalpha/Beta around the phosphorylation site of S180/181.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HepG2 cells using Phospho-IKKalpha/Beta (S180/181) Polyclonal Antibody.