

Phospho-JAK2 (Y1007) Antibody



PACO00227

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

JAK2

Uniprot

O60674

Synonyms:

JAK2; Tyrosine-protein kinase JAK2; Janus kinase 2; JAK-2


Immunogen:

Synthesized peptide derived from human JAK2 around the phosphorylation site of Y1007.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of NIH-3T3 cells using Phospho-JAK2 (Y1007) Polyclonal Antibody.