

Phospho-CDKN1B (S10) Antibody



PACO00248

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

CDKN1B

Uniprot

P46527

Synonyms:

CDKN1B; KIP1; Cyclin-dependent kinase inhibitor 1B; Cyclin-dependent kinase inhibitor p27; p27Kip1


Immunogen:

Synthesized peptide derived from human p27 around the phosphorylation site of S10.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using Phospho-p27 (S10) Polyclonal Antibody.