

Phospho-SHB (Y246) Antibody



PACO00266

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

SHB

Uniprot

Q15464

Synonyms:

SHB; SH2 domain-containing adapter protein B


Immunogen:

Synthesized peptide derived from human Shb around the phosphorylation site of Y246.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293 cells using Phospho-Shb (Y246) Polyclonal Antibody.