

Phospho-SP1 (T739) Antibody



PACO00267

Product Information

Size:

50ug

Protein Background:**Reactivity:**

Human, Mouse, Rat, Monkey

Gene ID:

SP1

Source:

Rabbit

Uniprot

P08047

Isotype:

IgG

Synonyms:

SP1; TSFP1; Transcription factor Sp1

Applications:

ELISA, WB, IHC

Immunogen:

Synthesized peptide derived from human Sp1 around the phosphorylation site of T739.


Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COS7 cells using Phospho-Sp1 (T739) Polyclonal Antibody.