

Phospho-STAT4 (Y693) Antibody



PACO00268

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IP

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300, IP:2-5ug

Protein Background:**Gene ID:**

STAT4

Uniprot

Q14765

Synonyms:

STAT4; Signal transducer and activator of transcription 4


Immunogen:

Synthesized peptide derived from human Stat4 around the phosphorylation site of Y693.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COS7 cells using Phospho-Stat4 (Y693) Polyclonal Antibody.