

Phospho-STAT5A/STAT5B (S726/731) Antibody



PACO00310

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

STAT5A/STAT5B

Uniprot

P42229/P51692

Synonyms:

STAT5A; STAT5; Signal transducer and activator of transcription 5A; STAT5B; Signal transducer and activator of transcription 5B


Immunogen:

Synthesized peptide derived from human Stat5 around the phosphorylation site of S726/731.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using Phospho-Stat5 (S726/731) Polyclonal Antibody.