

# Phospho-EIF2AK3 (T981) Antibody



PACO00329

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

IHC:1:100-1:300

**Protein Background:****Gene ID:**

EIF2AK3

**Uniprot**

Q9NZJ5

**Synonyms:**

EIF2AK3; PEK; PERK; Eukaryotic translation initiation factor 2-alpha kinase 3; PRKR-like endoplasmic reticulum kinase; Pancreatic eIF2-alpha kinase; HsPEK

**Immunogen:**


Synthesized peptide derived from human PERK around the phosphorylation site of T981.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of 453 AD293 cells using Phospho-PERK (T981) Polyclonal Antibody.