

Phospho-SNCA (Y125) Antibody



PACO00330

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

SNCA

Uniprot

P37840

Synonyms:

SNCA; NACP; PARK1; Alpha-synuclein; Non-A beta component of AD amyloid; Non-A4 component of amyloid precursor; NACP


Immunogen:

Synthesized peptide derived from human Synuclein-alpha around the phosphorylation site of Y125.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of SH-SY5Y 293 SK-HEP-1 cells using Phospho-Synuclein-alpha (Y125) Polyclonal Antibody.