

YWHAE Antibody



PACO00340

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

YWHAE

Uniprot

P62258

Synonyms:

YWHAE; 14-3-3 protein epsilon; 14-3-3E


Immunogen:

Synthesized peptide derived from the C-terminal region of human 14-3-3 epsilon.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using 14-3-3 ε Polyclonal Antibody.