

TP53BP1 Antibody



PACO00348

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:

Gene ID:

TP53BP1

Uniprot

Q12888

Synonyms:

TP53BP1; Tumor suppressor p53-binding protein 1; 53BP1; p53-binding protein 1; p53BP1

Immunogen:

Synthesized peptide derived from human 53BP1 around the non-phosphorylation site of S6.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of A549 cells using 53BP1 Polyclonal Antibody.

