

AARSD1 Antibody



PACO00350

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

AARSD1

Uniprot

Q9BTE6

Synonyms:

AARSD1; Alanyl-tRNA editing protein Aarsd1; Alanyl-tRNA synthetase domain-containing protein 1

Immunogen:

Synthesized peptide derived from the Internal region of human AARSD1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HepG2 cells using AARSD1 Polyclonal Antibody.