

# ABCC12 Antibody



PACO00354

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

ABCC12

**Uniprot**

Q96J65

**Synonyms:**

ABCC12; MRP9; Multidrug resistance-associated protein 9; ATP-binding cassette sub-family C member 12

**Immunogen:**

Synthesized peptide derived from the Internal region of human ABCC12.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of HuvEc cells using ABCC12 Polyclonal Antibody.

