

# ABHD11 Antibody



PACO00355

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

ABHD11

**Uniprot**

Q8NFV4

**Synonyms:**

ABHD11; WBSCR21; PP1226; Abhydrolase domain-containing protein 11; Williams-Beuren syndrome chromosomal region 21 protein

**Immunogen:**

Synthesized peptide derived from the Internal region of human ABHD11.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of HuvEc cells using ABHD11 Polyclonal Antibody.

