

EPHX4 Antibody



PACO00359

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

EPHX4

Uniprot

Q8IU55

Synonyms:

EPHX4; ABHD7; EPHXRP; Epoxide hydrolase 4; Abhydrolase domain-containing protein 7; Epoxide hydrolase-related protein

Immunogen:

Synthesized peptide derived from the Internal region of human ABHD7.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HT29 cells using ABHD7 Polyclonal Antibody.

