

ADCK1 Antibody



PACO00370

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

ADCK1

Uniprot

Q86TW2

Synonyms:

ADCK1; Uncharacterized aarF domain-containing protein kinase 1


Immunogen:

Synthesized peptide derived from the Internal region of human ADCK1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of RAW264.7 cells using ADCK1 Polyclonal Antibody.