

ADCK5 Antibody



PACO00372

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

ADCK5

Uniprot

Q3MIX3

Synonyms:

ADCK5; Uncharacterized aarF domain-containing protein kinase 5


Immunogen:

Synthesized peptide derived from the Internal region of human ADCK5.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HuvEc cells using ADCK5 Polyclonal Antibody.