

MLLT10 Antibody



PACO00374

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

MLLT10

Uniprot

P55197

Synonyms:

MLLT10; AF10; Protein AF-10; ALL1-fused gene from chromosome 10 protein


Immunogen:

Synthesized peptide derived from the Internal region of human AF-10.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using AF-10 Polyclonal Antibody.