

ALMS1 Antibody



PACO00388

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

ALMS1

Uniprot

Q8TCU4

Synonyms:

ALMS1; KIAA0328; Alstrom syndrome protein 1


Immunogen:

Synthesized peptide derived from the Internal region of human ALMS1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of various cells using ALMS1 Polyclonal Antibody.