

ANKRD30A Antibody



PACO00394

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

ANKRD30A

Uniprot

Q9BXX3

Synonyms:

ANKRD30A; Ankyrin repeat domain-containing protein 30A; Serologically defined breast cancer antigen NY-BR-1

Immunogen:

Synthesized peptide derived from the C-terminal region of human ANKRD30A.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of COLO205 cells using ANKRD30A Polyclonal Antibody.

