

NUSAP1 Antibody



PACO00395

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

NUSAP1

Uniprot

Q9BXS6

Synonyms:

NUSAP1; ANKT; BM-037; Nucleolar and spindle-associated protein 1; NuSAP


Immunogen:

Synthesized peptide derived from the C-terminal region of human ANKT.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HT29 cells using ANKT Polyclonal Antibody.